

## MW lamp LMW-15 (312 nm 20041)

### Set up conditions

The lamp was warmed up for >10 min.

Black felt paper was behind the lamp to block all reflections.

The OO (Fiber Optics Integrating Sphere) FOIS-1 was on the edge of the black felt paper about 7 ¼ in. from the surface of the lamp.

Fiber optic P1000-2-UV/VIS cable used to connect the FOIS and the USB200 spectrometer.

### Scan setting


Wavelength range is from 270 to 400 nm.

Integ. Time = **37**

Average = 2

Boxcar = 2

Flash Delay = 14



Measured 4/6/02

By DEN.